



Post Surgery Liposuction Compression



Since the beginning, the know-how, experience, and creativity in the manufacture of compression garments and corsetry was what we seek to produce the best product can be achieved.

Accuracy in the details, tissue quality, technical performance in line with the growing surgical requirements.

Overplast successfully finalizes all types of surgeries: facelift, tummy tuck, liposuction, breast augmentation, otoelastic.

The sheaths ensure constant and uniform skin compression throughout the post-operative course, offering a comfortable fit for the patient.

The use of advanced fabrics gives the patient a feeling of well-being, protection, and lightness.

First choice for the excellence of Italian cosmetic surgery.

Modern medicine uses compression therapies for the following reasons:


- To minimize swelling, bruises, and the occurrence of keloid scars.

- To reduce the risk of infection.

- To shape the contours and stabilize the implants.

- To reduce pain.

- To speed up the healing process.







Benefits

Defined ductility and elasticity
Maintain a stable microclimate environment
Antibacterial and hypoallergenic effect.

Eliminate the occurrence of unpleasant odor
Microfibers take excess moisture away from the body
Light and pleasant to the touch, making long-term wear easier.

Distinctive and soft texture constantly “massages” scars and thus speeds their healing. Latex-free. Flat seams do not put a strain on the operated parts.

It guarantees the required compression ranging between 20 ~ 30 mmHg.




BRA

Size	S	M	L	XL
Under Breast	69-82	79-92	89-102	102-120

SHEATHS

Size	40-42	44-46	48-50	52-54
Abdomen	67-69	77-89	87-99	99-118
Buttock	92-108	100-116	108-124	122-140
Haunch	54-64	61-71	68-78	77-97
Middle haunch	42-50	46-54	50-58	56-77
Under knee	30-37	34-41	38-45	43-65
Ankle	19-22	22-26	26-30	30-36

The measures in this table are in centimeters




LIPOSUCTION

The main properties are:

- The elasticity of compression and containment; the compression is carried out in all four directions to ensure containment in any position assumed by the patient.
- Resistance: the sheaths never lose their shape initial also after several washes; they are researched and designed to be h7 , 24 days on 7, then necessarily subject to continuous and frequent washings.
- High breathability and absolute non-allergenic; it does not alter the natural body moisture to evaporation processes: promotes maintaining the feeling of fresh and dry skin and minimizes the risk of irritation leading an itch and redness skin (issue encountered frequently in sheaths with traditional fabric). The hems at both ends of the sleeves, ensure the patient to enter but not contact chat with the outer fabric, ensuring high comfort and avoiding itching.

Lipoclyde sheaths present microfiber and other fabrics innovative grade to ensure maximum breathability, hypoallergenic, and while high comfort for the patient. The microfiber and a fabric composed of 58% from stretch nylon (family of polyamides synthetic) and 42% lycra (elastomeric fiber with unique properties of elongation and recovery of the original shape).

The sheath is invisible under clothes and provides a proper post-operative recovery.



Lipoclyde provides Panty Girdle Long and Short Ankle, foot including.



Panty short under the knee

Benefits: Easy donning with effective compression.
Closure available with hooks or velcro
Sizes: 38 - 40 - 42- 44 - 46 - 48 - 50 - 52 - 54 - 56
Composition: 77% polyamide, 23% elastane



All male and female garments with trouser legs

Panty short above knee

Benefits: Easy donning with effective compression.

Closure available with hooks or velcro

Sizes: 38 - 40 - 42- 44 - 46 - 48 - 50 - 52 - 54 - 56

Composition: 77% polyamide, 23% elastane



All male and female garments are available with closure in hooks or velcro.



Panty long ankle

Benefits: Easy donning with effective compression.

Sizes: 38 - 40 - 42- 44 - 46 - 48 - 50 - 52 - 54 - 56

Composition

LIPOSUCTION & ABDOMINOPLASTY

What is a post-surgical body?

The body post-surgery is an elasticized containment sheath used to reshape the body after an operation of liposculpture or liposuction. It allows the skin tissue to keep in touch with the muscular part and then get the silhouette desired.

There are allergy problems?

No, in fact, contact with the body of an elastic sheath gives an immediate sense of comfort and comfort. It does not irritate the skin because the inside the fabric is very smooth and soft.

What are the short anti-cellulite?

Anti-cellulite shorts are special elastic sleeves as within them employ silicone. With the normal movement of the legs, the silicone exerts a continuous massage which allows a reduction of cellulite localized on the thighs of 20%.



Liposuction & Abdominal models, allow you to remove the straps for greater daily comfort.

provides Panty Girdle Long and Short Ankle, foot including.



Mid-Thigh Girdle

Benefits: Easy donning with effective compression.

Sizes: 38 - 40 - 42- 44 - 46 - 48 - 50 - 52

Composition



Girdle short knee

Benefits: Easy donning with effective compression.

Sizes: 38 - 40 - 42 - 44 - 46 - 48 - 50 - 52

Composition





Girdle long knee


Benefits: Easy donning with effective compression.

Closure available with hooks or velcro

Sizes: 38 - 40 - 42 - 44 - 46 - 48 - 50 - 52

Composition: 77% polyamide, 23% elastane



 All male and female garments are available in three different colors: white, black, and skin.





Girdle long ankle

Benefits: Easy donning with effective compression.

Closure available with hooks or velcro

Sizes: 38 - 40 - 42- 44 - 46 - 48 - 50 - 52

Composition: 77% polyamide, 23% elastane

All girdles models are available with a removable or fixed shoulder strap to support the scar's medication

P PREGNANCY & ABDOMINOPLASTY

What is a post-surgical body?

The body post-surgery is an elasticized containment sheath used to reshape the belly after pregnancy or liposculpture to reduce the risk of stretch-marks and scar's keloids.

It allows the skin tissue to keep in touch with the muscular part and then get the silhouette desired.

There are allergy problems?

No, in fact, contact with the body of an elastic sheath gives an immediate sense of comfort and comfort. It does not irritate the skin because the inside the fabric is very smooth and soft.



Abdominal girdle

Benefits: Easy donning with effective compression with adjustable velcro closure

Sizes: S - M - L - XL - XXL

Composition





Abdominal band girdle

Benefits: Easy and fast donning with effective compression with modurable two or three bands

Sizes: S - M - L - XL - XXL

Composition: 25% polyamide, 10% cotton, 65% elastane

Abdominal long girdle

Benefits: Easy donning with effective compression with adjustable velcro closure

Sizes: S - M - L - XL - XXL

Composition: 25% polyamide, 10% cotton, 65% elastane





Culotte

Benefits: Easy donning with effective compression.

Sizes: 38 - 40 - 42- 44 - 46 - 48 - 50 - 52

Composition


B_{RA}

What are the post-surgical bras?

The post-surgical bras keep the breasts up, separated, and in the form that the surgeon wanted to impress. The innovative front connection through hooks allows you to adjust the voltage and then to decrease pain after augmentation. It also helps rest better in the hours of sleep.

How long you need to wear a post-surgical sleeve after the operation?

Normally we can say that the best results in a surgical operation are reached three months after surgery. We say so that at this time would be indicated always wear the sheath.



Overbra

Benefits: Versatile model with an elastic band to be positioned according to the post-surgical needs. Stabilizes the implant

Sizes: S - M - L - XL

Composition: 59% polyamide, 22% cotton, 19% elastane



Oversplit

Benefits: Stabilizes the prosthetic implant to remain free of the space between the breasts

Sizes: S - M - L - XL

Composition: 59% polyamide, 22% cotton, 19% elastane



Overshape

Benefits: Stabilizes the prosthetic implant. Comfortable closure of velcro straps

Sizes: S - M - L - XL

Composition: 59% polyamide, 22% cotton, 19% elastane





Implant support bra

Benefits: Stabilize the prosthesis ensuring a high compression of tissues, supported by the support of an underbust band

Sizes: S - M - L - XL

Composition: 77% polyamide, 23% elastane



Support asymmetric bra

Benefits: Indicated for breast asymmetries, supporting high-compression straps with Velcro closures comfortable

Sizes: S - M - L - XL

Composition: 77% polyamide, 23% elastane




Strapless

Benefits: Stabilizes the prosthetic implant with the aid of an adjustable elastic band. Model strapless

Sizes: S - M - L - XL

Composition: 77% polyamide, 23% elastane



A woman is shown from the waist up, wearing a black and white Double Stabilizer Band. The band consists of two horizontal bands, one black and one white, with a white band in the center. The background is a brick wall.

Double Stabilizer Band

Benefits: Two stabilizing bands with adjustable closure according to the post-surgical requirements

Sizes: S - M - L - XL

Composition: 5% polyamide, 5% cotton, 90% elastane



A woman is shown from the waist up, wearing a white Implant stabilizer band. The band is a single horizontal band with a white band in the center. The background is a brick wall.

Implant stabilizer band

Benefits: One band with adjustable closure according to the post-surgical requirements

Sizes: S - M - L - XL

Composition: 5% polyamide, 5% cotton, 90% elastane



- All bra is available in three different colors: white, black and skin.
The most of all bra is available with a shoulder strap.

MASKS

The innovative design and easy to use with an adjustable strap.

For a facelift, liposuction of double chin, and otoplasty. The wide lateral velcro straps mean the support can be altered to suit the type of intervention and the patient's shape.

The fabric will not run, this allows the physician to make an incision over the ear so the support can be maintained.





Otoelastic

Benefits: High compression with adjustable velcro closure

Size: Standard

Composition: 77% polyamide, 23% elastane



Face Lifting



Face lifting

Benefits: High compression with adjustable velcro closure

Size: Standard

Composition: 77% polyamide, 23% elastane



Mentelastic



Mentoelastic

Benefits: High compression with adjustable velcro closure

Size: Standard


Composition: 77% polyamide, 23% elastane



SPECIAL MODELS

Our laboratory can meet any surgical needs manufacturing a specific model with special features.

We can use a different type of texture according to surgeon demand.



Brassieres

Benefits: Full garment to meet the mastoplasty and arm liposuction's needs

Sizes: S - M - L - XL

Composition: 77% Polyamide, 23% elastane



Overarm

Benefits: It is indicated for patients who have had an arm lift. Excellent comfort and support

Sizes: S - M - L - XL

Composition: 77% Polyamide, 23% elastane

